

Aidaptus® Specifications

Needle safety	Automatic needle shielding with needle hidden before, during & after use	Needle insertion	Automatic
Primary container/ syringe type	<p>1mL small round and cut/cropped flange (ISO standard)</p> <ul style="list-style-type: none"> BD Hypak, BD Neopak, Gerresheimer, Nipro, Ompi <p>2.25 mL small round flange (ISO standard)</p> <ul style="list-style-type: none"> BD Hypak, BD Neopak, Ompi. 	Start and end of dose notification	End of dose click (>70Db at 1m)
Fill Volume	<ul style="list-style-type: none"> 1mL: 0.3 - 1.0mL * 2.25mL: 0.5 - 2.0mL* Adapts to a range of different drug fill volumes with no changed parts, with a self-adjusting plunger. <p><small>*includes air bubble and overfill</small></p>	Weight	28g (without syringe)
Needle type and gauge	27G STW, 29G TW	Dimensions	H = 162mm (cap on) W=18mm
Needle insertion depth (subcutaneous)	1mL = 6.5mm 2.25 mL = 6.3mm Syringe dependent	Component count	<p>1mL = 11 plastic components, 3 springs, 2 metal pressings</p> <p>2.25mL = 10 plastic components, 3 springs, 2 metal pressings</p>
Injection angle	90°		